

Package ‘PKbioanalysis’

January 14, 2026

Type Package

Title Pharmacokinetic Bioanalysis Experiments Design and Exploration

Version 0.4.0

Maintainer Omar Elashkar <omar.i.elashkar@gmail.com>

Description Automate pharmacokinetic/pharmacodynamic bioanalytical procedures based on best practices and regulatory recommendations.

The package impose regulatory constrains and sanity checking for common bioanalytical procedures.

Additionally, 'PKbioanalysis' provides a relational infrastructure for plate management and injection sequence.

License AGPL (>= 3)

Encoding UTF-8

RoxygenNote 7.3.3

Depends R (>= 4.3.0)

Imports dplyr, tidyselect, stringr (>= 1.5.1), ggplot2 (>= 3.5.1), ggforce (>= 0.4.1), tidyr, glue (>= 1.6.2), checkmate, shiny, DBI, duckdb (>= 1.0.0), bslib, bsicons, cli, DiagrammeR, DT, shinyWidgets, shinyjs, units, stats, shinyalert, htmltools, rlang, grDevices, utils, yaml, rhandsontable, methods, reticulate, RTMB, janitor, tibble, xml2, RaMS, data.tree, nloptr, forcats, ggiraph, gt, plotly, pracma, reactable, gtools, ellmer, htmlwidgets, jsonlite, nlme, pmxTools, scales, shinychat, uuid, sortable

Suggests rxode2, knitr, patchwork, rmarkdown, testthat (>= 3.0.0)

SystemRequirements python (>= 3.10)

URL <https://omarashkar.github.io/PKbioanalysis/>

BugReports <https://github.com/OmarAshkar/PKbioanalysis/issues>

Config/testthat/edition 3

Collate 'PKbioanalysis-package.R' 'calc.R' 'chrom_anomaly_plots.R' 'chrom_app.R' 'chrom_parsers.R' 'chrom_utils.R' 'generics.R' 'class.R' 'chromatogram.R' 'config.R' 'dil_map.R'

```
'estim_residuals.R' 'genAI.R' 'gen_studydesign.R'
'generate_test_set.R' 'injec_list.R' 'linearitycheck.R'
'method_file.R' 'peak_integrate.R' 'plate.R' 'plate_expr.R'
'plates_class.R' 'process_chroms.R' 'quant_app.R'
'quant_misc.R' 'quant_object.R' 'quant_parsers.R'
'samples_profiles.R' 'study_app.R' 'suitability_utils.R'
'utils.R' 'writers.R' 'zzz.R'
```

NeedsCompilation no

Author Omar Elashkar [aut, cre] (ORCID:

<<https://orcid.org/0000-0002-5505-778X>>)

Repository CRAN

Date/Publication 2026-01-14 20:40:07 UTC

Contents

add_blank	3
add_cs_curve	4
add_DB	5
add_DQC	5
add_QC	6
add_samples	7
add_samples_db	8
add_samples_db2	8
add_suitability	9
area_report.PeakRes	10
calc_var_summary	10
check_chrom_cmpds	11
chrom_app	12
combine_injec_lists	12
combine_plates	13
config_suitability	13
create_new_study	14
cv	14
download_sample_list	15
estim_dil_limit	15
estim_lloq	16
export_integration	16
export_run	17
extract_peak_bounds	17
fill_scheme	18
filter_chrom	18
fit_var	19
formatted_print	20
generate_96	20
get_compound_ID	21
get_sample_ID	21

get_sample_names	22
has_default_RT	22
install_py_dep	23
integrate	23
is_integrated	24
is_smoothed	24
length,MultiPlate-method	25
make_calibration_study	25
make_metabolic_study	26
plate_metadata	27
plate_tree	27
plot.PlateObj	28
plot_chrom	29
plot_peak_areas.PeakRes	30
plot_RT.ChromRes	31
plot_RT.PeakRes	31
plot_var_pattern	32
precision_per_vial	32
prefilter_precision_data	33
quant_app	34
read_chrom	34
read_experiment_results	35
register_plate	35
response_to_conc	36
reverse_predict	36
run_summary	37
smooth_chrom	37
study_app	38
update_RT	38
write_injec_seq	39
[[,MultiPlate-method	39
Index	41

add_blank	<i>Add blank to the plate Can be either double blank (DB), CS0IS+ or CS+IS0</i>
-----------	---

Description

Add blank to the plate Can be either double blank (DB), CS0IS+ or CS+IS0

Usage

```
add_blank(plate, IS = TRUE, analyte = FALSE, analytical = FALSE, group = NA)
```

Arguments

plate	PlateObj object
IS	logical. If TRUE, add IS to the well.
analyte	logical. If TRUE, add analyte to the well.
analytical	logical. If FALSE, the blank is analytical, if TRUE it is bioanalytical.
group	A string for bioanalytical group.

Value

PlateObj

add_cs_curve	<i>Add calibration curve to the plate</i>
--------------	---

Description

Add calibration curve to the plate

Usage

```
add_cs_curve(plate, plate_std, rep = 1, group = NA)
```

Arguments

plate	PlateObj
plate_std	character
rep	numeric. Number of technical replicates. Default is 1.
group	A string for bioanalytical group.

Value

PlateObj

Examples

```
plate <- generate_96() |>
  add_cs_curve(c(1, 3, 5, 10, 50, 100, 200))
plot(plate)
```

add_DB	<i>Add double blank (DB) to a plate</i>
--------	---

Description

Add double blank (DB) to a plate

Usage

```
add_DB(plate, analytical = FALSE, group = NA)
```

Arguments

plate	PlateObj object
analytical	logical. If TRUE, the blank is bioanalytical, if FALSE it is analytical.
group	A string for bioanalytical group.

Value

PlateObj

Examples

```
plate <- generate_96() |>  
add_DB()
```

add_DQC	<i>Add dilution quality control (DQC) to the plate</i>
---------	--

Description

Add dilution quality control (DQC) to the plate

Usage

```
add_DQC(plate, conc, fac, rep = 5, group = NA)
```

Arguments

plate	PlateObj object
conc	numeric. Concentration of the DQC well.
fac	numeric. Factor of the DQC well.
rep	numeric. Number of replicates. Default is 5.
group	A string for bioanalytical group. The current implementation does not check ULOQ or LLOQ boundaries.

`add_QC`*Add quality control samples to the plate*

Description

A function to add QCs to plate. This function assumes adherence to ICH guideline M10 on bioanalytical method validation and study sample analysis Geneva, Switzerland (2022). If you are not following this guideline, you can set 'reg = TRUE' to ignore the restrictions.

Usage

```
add_QC(  
  plate,  
  lqc_conc,  
  mqc_conc,  
  hqc_conc,  
  extra = NULL,  
  n_qc = 3,  
  qc_serial = TRUE,  
  reg = TRUE,  
  group = NA  
)
```

Arguments

<code>plate</code>	PlateObj object
<code>lqc_conc</code>	low quality control concentration
<code>mqc_conc</code>	medium quality control concentration
<code>hqc_conc</code>	high quality control concentration
<code>extra</code>	numeric vector of extra QC concentrations.
<code>n_qc</code>	number of QC sets. Default is 3
<code>qc_serial</code>	logical. If TRUE, QCs are placed serially
<code>reg</code>	logical. Indicates if restrictions should not be applied to the QC samples. Default is TRUE
<code>group</code>	A string for bioanalytical group.

Value

PlateObj

add_samples

Add samples to plate with pharmacokinetic attributes

Description

Add samples to plate with pharmacokinetic attributes

Usage

```
add_samples(plate, samples, prefix = NA, dil = NA, group = NA, rep = 1)
```

Arguments

plate	PlateObj
samples	A vector representing samples names. Must be unique.
prefix	A prefix to be added before samples names. Default is "Sub".
dil	A vector representing samples' dilution factor. Must be same length as samples.
group	A vector representing samples' bioanalytical group. Must be same length as samples.
rep	Number of technical replicates for each combination. Default is 1.

Details

final name will be of form. Prefix-SampleName-Time-Concentration-Factor samples must be a unique vector and did not exist in the plate before. Time is either a vector or a single value. If it is a vector, it will be repeated for each sample. Conc, dil, factor and dose are either a vector or a single value. If it is a vector, it must be the corresponding length of samples.

Allowed routes are "IV", "IM", "IP", "SC", "PO", "INH" which are short for Intravenous, Intramuscular, Intraperitoneal, Subcutaneous, Per Os (oral), Inhalation.

Factor is an arbitrary factor used in the design like food vs fasted, healthy vs diseased, positive genotype ... etc.

Value

PlateObj

Examples

```
plate <- generate_96() |>
  add_samples(paste0("T", 1:12))
```

add_samples_db	<i>Add samples from the sample log to the plate</i>
----------------	---

Description

Add samples from the sample log to the plate

Usage

```
add_samples_db(plate, logIds, dil = 1, namestyle = 1, group = NA)
```

Arguments

plate	PlateObj
logIds	A vector of log IDs from the sample log.
dil	A vector with length corresponding number of logIds. See details.
namestyle	A numeric value indicating the naming style. 1 for long names, 2 for short names.
group	A string for bioanalytical group.

Details

This function will retrieve sample information from the sample log database using the provided log IDs. It constructs sample names based on the specified naming style and adds them to the plate. The 'dil' parameter allows specifying dilution factors for each sample, which will be appended to the sample names. If a single dilution factor is provided, it will be applied to all samples.

Value

PlateObj

add_samples_db2	<i>Add samples from the sample log to the plate with multiplication</i>
-----------------	---

Description

Add samples from the sample log to the plate with multiplication

Usage

```
add_samples_db2(plate, logIds, dil = c(1, 1), namestyle = 1, group = NA)
```


Arguments

plate	PlateObj
logIds	A vector of log IDs from the sample log.
dil	A vector with length corresponding number of repeats. See details.
namestyle	A numeric value indicating the naming style. 1 for long names, 2 for short names.
group	A string for bioanalytical group.

Details

This function is wrapper around 'add_samples_db()' that allows for quick replication of samples by dilution factor vector. For instance, if `dil = c(1,10)`, the samples will be repeated twice with one fold and 10 fold dilution factor each time.

add_suitability	<i>Add suitability sample to the plate</i>
-----------------	--

Description

Add suitability sample to the plate

Usage

```
add_suitability(plate, conc, label = "suitability", group = NA)
```

Arguments

plate	PlateObj object.
conc	numeric. Concentration of the suitability well.
label	character. Label for the suitability well. Default is "suitability".
group	A string for bioanalytical group.

Value

PlateObj

area_report.PeakRes *gt table of areas*

Description

gt table of areas

Usage

```
area_report.PeakRes(
  peaks_res,
  normalize = TRUE,
  blanks = TRUE,
  analytes = TRUE,
  standards = TRUE,
  QCs = TRUE,
  compounds = NULL
)
```

Arguments

peaks_res	PeakRes object
normalize	logical. If TRUE, normalize the peak area by the IS area.
blanks	logical. If TRUE, include blanks
analytes	logical. If TRUE, include analytes
standards	logical. If TRUE, include standards
QCs	logical. If TRUE, include QCs
compounds	numeric vector of compound numbers to include. If NULL, include all compounds

calc_var_summary	<i>Calculate Summary Statistics for Each Concentration Level For Either Concentration, Area, or Area Ratio</i>
------------------	--

Description

Calculate Summary Statistics for Each Concentration Level For Either Concentration, Area, or Area Ratio

Usage

```
calc_var_summary(  
  df,  
  col = "conc",  
  acc_cutoff = 0.2,  
  dev_cutoff = 0.2,  
  type = "QC"  
)
```

Arguments

df	Data frame with columns: stdconc (standardized concentration), conc (concentration), area (peak area), area_ratio (area ratio)
col	Column to calculate summary for ("conc", "area", or "area_ratio")
acc_cutoff	Accuracy threshold (default is 20%) for concentration vs standard concentration
dev_cutoff	Deviation threshold (default is 20%) for concentration vs standard concentration
type	Type of samples to include ("Standard", "QC", "DQC")

Author(s)

Omar I. Elashkar

check_chrom_cmpds	<i>Check Matching of Compound and Transitions in chrom_res and method database</i>
-------------------	--

Description

Check Matching of Compound and Transitions in chrom_res and method database

Usage

```
check_chrom_cmpds(chrom_res, method_id)
```

Arguments

chrom_res	ChromRes object
method_id	Method ID in the method database This is important to give no error before merging quantification results to ensure consistency.

Value

TRUE if all compounds and transitions match, otherwise FALSE

chrom_app	<i>chrom_apps</i>
-----------	-------------------

Description

This function creates a shiny app for peak integration.

Usage

```
chrom_app()
```

combine_injec_lists	<i>Create Sample List with rigorous design</i>
---------------------	--

Description

Create Sample List with rigorous design

Usage

```
combine_injec_lists(  
  sample_lists,  
  n_equi = 10,  
  equi_pos,  
  equi_prefix = Sys.Date(),  
  equi_suffix = "equi",  
  equi_injec_vol = 0.5  
)
```

Arguments

sample_lists	a list of sample lists
n_equi	number of equilibration injections
equi_pos	position of equilibration injections. For format check details
equi_prefix	prefix for equilibration injections
equi_suffix	suffix for equilibration injections
equi_injec_vol	volume of equilibration injection

Details

The equi_pos format will be Row:Column format. E.g: "A,1"

Value

InjecListObj object

combine_plates	<i>Combine plates in MultiPlate object</i>
----------------	--

Description

Combine plates in MultiPlate object

Usage

```
combine_plates(plates)
```

Arguments

plates list of PlateObj objects

Value

MultiPlate object

config_suitability	<i>Configure suitability runs</i>
--------------------	-----------------------------------

Description

Configure suitability runs by specifying vial position and range of runs to include.

Usage

```
config_suitability(quantres, vial_pos, start = NULL, end = NULL)
```

Arguments

quantres QuantRes object

vial_pos Vial position to use for suitability (e.g., "2:H,9")

start Start position (1-based index) of runs to include. If NULL, starts from the first run.

end End position (1-based index) of runs to include. If NULL, ends at the last run.

Value

Updated QuantRes object with suitability configuration.

create_new_study	Create a new study in the database
------------------	------------------------------------

Description

Create a new study in the database

Usage

```
create_new_study(df)
```

Arguments

df	A data frame with one row containing study details: type, title, description, pkstudy (logical), subject_type
----	---

Value

A data frame with the created study details including generated id, start_date, and status

cv	Calculate Coefficient of variation
----	------------------------------------

Description

Calculate Coefficient of variation

Usage

```
cv(x, percent = TRUE)
```

Arguments

x	vector
percent	To return the value as percentage

Details

A simple calculation of the coefficient of variation (CV) is done as the standard deviation divided by the mean. By default, the result is in percentage.

Value

numeric

download_sample_list	<i>Download sample list from database to local spreadsheet with vendor specific format</i>
----------------------	--

Description

Download sample list from database to local spreadsheet with vendor specific format

Usage

```
download_sample_list(sample_list, vendor)
```

Arguments

sample_list	dataframe of sample list either from db or from write_injec_seq
vendor	currently only 'masslynx', 'masshunter' and 'analyst' are supported

Details

For all current vendors, the exported format will be in csv format, compatible with the respective software.

Value

dataframe

estim_dil_limit	<i>Estimate Dilution Limit Based on Additive and Proportional Errors and LLOQ</i>
-----------------	---

Description

Estimate Dilution Limit Based on Additive and Proportional Errors and LLOQ

Usage

```
estim_dil_limit(add_err, prop_err, lloq)
```

Arguments

add_err	Additive error (constant)
prop_err	Proportional error (CV)
lloq	Lower limit of quantification

Author(s)

Omar I. Elashkar

Examples

```
estim_dil_limit(add_err=0.1, prop_err=0.1, lloq=1)
estim_dil_limit(add_err=1, prop_err=0.1, lloq=55)
```

estim_lloq

*Estimate LLOQ From Existing Additive and Proportional errors***Description**

Estimate LLOQ From Existing Additive and Proportional errors

Usage

```
estim_lloq(add_err = 0.04, prop_err = 0.05, cv_lloq = 0.2, cv_lqc = 0.15)
```

Arguments

add_err	Additive error (constant)
prop_err	Proportional error (CV)
cv_lloq	Maximum coefficient of variation at LLOQ
cv_lqc	Maximum coefficient of variation at LQC

A method to estimate LLOQ from existing additive and proportional errors. The function does inequality constrained optimization to find the LLOQ.

Author(s)

Omar I. Elashkar

export_integration

*Export Expected RT***Description**

Export expected RT values for each peak in the chromatogram.

Usage

```
export_integration(chrom_res, path)
```

Arguments

chrom_res	ChromRes object
path	Path to save the file

export_run	<i>Export run</i>
------------	-------------------

Description

Export run

Usage

export_run(peaks_res, path)

Arguments

peaks_res	PeakRes object
path	path to save csv

extract_peak_bounds	<i>Extract Peak Boundaries</i>
---------------------	--------------------------------

Description

Extract peak boundaries for a given compound ID

Usage

extract_peak_bounds(chrom_res, compound_id)

Arguments

chrom_res	ChromRes object
compound_id	Compound ID The function automatically prioritizes observed peak boundaries (manual integration) over expected ones. If observed boundaries are not available, it falls back to expected boundaries.

Value

Dataframe with compound_id, min, max

fill_scheme	<i>Filling orientation of the plate</i>
-------------	---

Description

This function sets the filling scheme of the plate. The filling scheme is used to determine the order in which the samples are filled in the plate. The default filling scheme is horizontal, which means that the samples are filled from left to right and top to bottom. The vertical filling scheme means that the samples are filled from top to bottom and left to right.

Usage

```
fill_scheme(  
  plate,  
  fill = "h",  
  tbound = "A",  
  bbound = "H",  
  lbound = 1,  
  rbound = 12  
)
```

Arguments

plate	PlateObj
fill	character. Filling scheme. Either "h" for horizontal, "v" for vertical.
tbound	character. Top bound of the filling scheme. Default is "A"
bbound	character. Bottom bound of the filling scheme. Default is "H"
lbound	numeric. Left bound of the filling scheme. Default is 1
rbound	numeric. Right bound of the filling scheme. Default is 12

Value

PlateObj

filter_chrom	<i>title Filter Chromatogram Peaks</i>
--------------	--

Description

This function filters chromatogram peaks based on transition ID and sample ID.

Usage

```
filter_chrom(  
  chrom_res,  
  transitions_ids = NULL,  
  samples_id = NULL,  
  cmpd_ids = NULL  
)
```

Arguments

chrom_res	ChromRes object
transitions_ids	Vector of transition IDs to filter. If NULL, all transitions are returned.
samples_id	Sample ID to filter.
cmpd_ids	Compound ID to filter. It must be numeric. If NULL, all compounds are returned.

fit_var	<i>Estimate Additive and proportional errors from calibration data</i>
---------	--

Description

Estimate Additive and proportional errors from calibration data

Usage

```
fit_var(  
  data,  
  level = 0.95,  
  method = "nlminb",  
  bootstrap = FALSE,  
  n_boot = 1000  
)
```

Arguments

data	Data frame with columns: conc (concentration), stdconc (standardized concentration, e.g. conc/LLOQ)
level	Confidence level for the CI (default is 0.95)
method	Optimization method (default is "nlminb")
bootstrap	Logical indicating whether to perform bootstrap (default is TRUE)
n_boot	Number of bootstrap samples (default is 1000)

Author(s)

Omar I. Elashkar

formatted_print	<i>Format and print the results of fit_var</i>
-----------------	--

Description

Format and print the results of fit_var

Usage

```
formatted_print(x, digits = 3)
```

Arguments

- x Data frame with results
- digits Number of digits to display

Author(s)

Omar I. Elashkar

generate_96	<i>Generate 96 well plate</i>
-------------	-------------------------------

Description

Generate 96 well plate

Usage

```
generate_96(descr = "", start_row = "A", start_col = 1)
```

Arguments

- descr plate description.
- start_row A letter corresponding to empty rows in a 96 well plate. Default is A.
- start_col A number indicating a column number to start with, given the start row. Default is 1.
Generate a typical 96 well plate. User need to specify the empty rows which a going to be used across the experiment.

Value

PlateObj

Examples

```
plate <- generate_96()
plot(plate)

plate <- generate_96("calibration", start_row = "C", start_col = 11)
plot(plate)
```

get_compound_ID	<i>Find Compound ID from compound Name</i>
-----------------	--

Description

This function returns the compound ID

Usage

```
get_compound_ID(chrom_res, compound_name)
```

Arguments

chrom_res	ChromRes object
compound_name	Compound Name

get_sample_ID	<i>Find Sample ID from sample Name</i>
---------------	--

Description

This function returns the sample ID

Usage

```
get_sample_ID(chrom_res, sample_name)
```

Arguments

chrom_res	ChromRes object
sample_name	Sample Name

get_sample_names	<i>Find sample names for all samples</i>
------------------	--

Description

Find sample names for all samples

Usage

```
get_sample_names(chrom_res)
```

Arguments

chrom_res	ChromRes object
-----------	-----------------

Value

data.frame with sample and sample_id

has_default_RT	<i>check if default expected RT is set for a compound</i>
----------------	---

Description

check if default expected RT is set for a compound

Usage

```
has_default_RT(chrom_res, compound_id)
```

Arguments

chrom_res	ChromRes object
compound_id	Compound ID

install_py_dep	<i>Install Python dependencies for PKbioanalysis</i>
----------------	--

Description

Install Python dependencies for PKbioanalysis

Usage

```
install_py_dep(..., envname = "PKbioanalysis")
```

Arguments

...	Additional arguments passed to reticulate::py_install
envname	Name of the virtual environment to create/use. Default is "PKbioanalysis"

Value

None

integrate	<i>integrate Peak with trapzoid method given start and end</i>
-----------	--

Description

integrate Peak with trapzoid method given start and end

Usage

```
integrate(chrom_res, compound_id, samples_ids, smoothed = TRUE)
```

Arguments

chrom_res	ChromRes object. Must have observed RT values
compound_id	Compound ID
samples_ids	Sample ID. If NULL, all samples will be used
smoothed	Logical. If TRUE, use smoothed chromatogram. Default is TRUE

is_integrated	<i>Check if peak was integrated for a specific compound</i>
---------------	---

Description

Check if peak was integrated for a specific compound

Usage

```
is_integrated(chrom_res, compound_id, sample_id = NULL)
```

```
## S4 method for signature 'ChromRes'  
is_integrated(chrom_res, compound_id, sample_id = NULL)
```

```
## S4 method for signature 'ChromResBase'  
is_integrated(chrom_res, compound_id, sample_id = NULL)
```

Arguments

chrom_res	ChromRes or ChromResBase object
compound_id	Compound ID
sample_id	Sample ID. If NULL, all samples are checked

Value

logical

Examples

```
## Not run:  
lapply(1:10, \(x) is_integrated(chrom_res, sample_id = 1, compound_id = 1))  
  
## End(Not run)
```

is_smoothed	<i>Return an indicator if the chromatogram is smoothed</i>
-------------	--

Description

Return an indicator if the chromatogram is smoothed

Usage

```
is_smoothed(chrom_res)
```

Arguments

chrom_res	ChromRes object
-----------	-----------------

length,MultiPlate-method

Length method for MultiPlate

Description

Length method for MultiPlate

Usage

```
## S4 method for signature 'MultiPlate'
length(x)
```

Arguments

x MultiPlate object

Value

number of plates

make_calibration_study

Create a calibration study with calibration standards and QCs

Description

Create a calibration study with calibration standards and QCs

Usage

```
make_calibration_study(
  plate,
  plate_std,
  lqc_conc = NULL,
  mqc_conc = NULL,
  hqc_conc = NULL,
  n_qc = NULL,
  qc_serial = FALSE,
  n_CS0IS0 = 1,
  n_CS0IS1 = 2,
  n_CS1IS0 = 1,
  group = NA
)
```

Arguments

plate	PlateObj object
plate_std	vector of calibration standards
lqc_conc	LQC concentration
mqc_conc	MQC concentration
hqc_conc	HQC concentration
n_qc	number of QC sets
qc_serial	logical. If TRUE, QCs are placed serially
n_CS0IS0	number of CS0IS0 (double) blanks
n_CS0IS1	number of CS0IS1 blanks
n_CS1IS0	number of CS1IS0 blanks
group	A string for bioanalytical group.

Value

PlateObj

make_metabolic_study *Create a metabolic study layout*

Description

Create a metabolic study layout

Usage

```
make_metabolic_study(
  study = "Metabolic Study",
  cmpds,
  time_points = c(0, 5, 10, 15, 30, 45, 60, 75, 90, 120),
  dose = NA,
  n_NAD = 3,
  n_noNAD = 2
)
```

Arguments

study	study name
cmpds	vector of compounds, including any standards
time_points	vector of time points
dose	dose amount. Default is NA
n_NAD	number of NAD positive samples. Default is 3
n_noNAD	number of NAD negative samples. Default is 2

Details

Note that this function does not require plate object. It will create a plate object automatically and return MultiPlate object

Value

MultiPlate object

plate_metadata	<i>Set plate description</i>
----------------	------------------------------

Description

Set plate description

Usage

```
plate_metadata(plate, descr)
```

Arguments

plate	PlateObj
descr	character. Description of the plate

Value

PlateObj

plate_tree	<i>Plot the design of the plate</i>
------------	-------------------------------------

Description

Plot the design of the plate

Usage

```
plate_tree(plate, plot = TRUE)
```

Arguments

plate	PlateObj object
plot	logical. If TRUE, plot the tree

Value

data.tree Node object or DiagrammeR object plot_tree will focus only on bioanalytical vial types, namely blanks, analytes, standards, QCs. The tree order will be plate_id, then group, then vial type, then entity, then number of technical replicates.

plot.PlateObj	<i>Plotting 96 well plate</i>
---------------	-------------------------------

Description

Plotting 96 well plate

Usage

```
## S3 method for class 'PlateObj'
plot(
  x,
  color = "conc",
  Instrument = "",
  caption = "",
  label_size = 1,
  transform_dil = FALSE,
  watermark = "auto",
  layoutOverlay = FALSE,
  path = NULL,
  ...
)
```

Arguments

x	PlateObj
color	character. Coloring variable. Choices: "conc", "group", "dil", "study", "time", "factor", "samples", "arm", "sex", "dose", "route", "matrix". Default is "conc"
Instrument	A string placed at subtitle
caption	A string place at plate caption
label_size	numeric. Size of the label. Default is 15
transform_dil	logical. If TRUE, transform the dilution factor to the label
watermark	character. If "auto", a watermark is added to the plot. If "none", no watermark is added. Default is "auto"
layoutOverlay	logical. If TRUE, overlay the plot layout. Default is FALSE
path	If not null, must be a path to save plate image
...	additional arguments passed to ggplot2::ggsave

Value

ggplot object

Examples

```
plate <- generate_96("new_plate", "C", 11) |>
  add_blank(IS = FALSE, analyte = FALSE) |>
  add_blank(IS = TRUE, analyte = FALSE) |>
  add_samples(c(
    "RD_per1", "RD_in1", "RD_T30", "RD_T60", "RD_T90", "RD_per2", "RD_in2",
    "EE_in0", "EE_T30", "EE_in30", "EE_T60", "EE_in60", "EE_T90", "EE_in90"
  ))
plot(plate)
```

plot_chrom

Plot Chromatogram per Sample for Selected transitions

Description

This function plots chromatograms for selected transitions per sample.

Usage

```
plot_chrom(
  chrom_res,
  ncol = 2,
  transitions_ids = NULL,
  sample_id,
  integrated = FALSE,
  show_RT = FALSE,
  smoothed = FALSE
)
```

Arguments

chrom_res	ChromRes object
ncol	Number of columns for facet_wrap. If 0, the chromatograms are overlayed in a single plot.
transitions_ids	Vector of transition IDs to plot. If NULL, all transitions are plotted.
sample_id	Sample ID to plot.
integrated	Boolean to show integrated area overlayed
show_RT	Boolean to show RT values
smoothed	Boolean to show smoothed chromatogram

Examples

```
## Not run:
path <- system.file("extdata", "waters_raw_ex", package="PKbioanalysis")
main <- read_chrom(path, method = 1)
plot_chrom(main, ncol = 2, transitions_ids = c(18,19,20), sample_id = 3)
plot_chrom(main, ncol = 3, transitions_ids = c(18,19,20), sample_id = 3)
plot_chrom(main, ncol = 1, transitions_ids = c(18,19,20), sample_id = 3)
plot_chrom(main, ncol = NULL, transitions_ids = c(18,19,20), sample_id = 3)

## End(Not run)
```

plot_peak_areas.PeakRes

Plot peak areas

Description

Plot peak areas

Usage

```
plot_peak_areas.PeakRes(
  peaks_res,
  normalize = TRUE,
  blanks = TRUE,
  compounds = NULL,
  analytes = TRUE,
  standards = TRUE,
  QCs = TRUE,
  type = "bar"
)
```

Arguments

peaks_res	PeakRes object
normalize	logical. If TRUE, normalize the peak area by the IS area.
blanks	logical. If TRUE, plot blanks
compounds	numeric vector of compound numbers to include. If NULL, include all compounds
analytes	logical. If TRUE, plot analytes
standards	logical. If TRUE, plot standards
QCs	logical. If TRUE, plot QCs
type	character. Either "bar" or "line"

Value

ggplot2 object

plot_RT.ChromRes	<i>Plotting RT intervals of chromatogram</i>
------------------	--

Description

Plotting RT intervals of chromatogram

Usage

```
plot_RT.ChromRes(chrom_res)
```

Arguments

chrom_res	ChromRes object
-----------	-----------------

plot_RT.PeakRes	<i>Plot RT</i>
-----------------	----------------

Description

Plot RT

Usage

```
plot_RT.PeakRes(  
  peaks_res,  
  normalize = TRUE,  
  blanks = TRUE,  
  analytes = TRUE,  
  standards = TRUE,  
  QCs = TRUE,  
  facet = FALSE,  
  compounds = NULL  
)
```

Arguments

peaks_res	PeakRes object
normalize	logical. If TRUE, normalize the peak area by the IS area.
blanks	logical. If TRUE, plot blanks
analytes	logical. If TRUE, plot analytes
standards	logical. If TRUE, plot standards
QCs	logical. If TRUE, plot QCs
facet	logical. If TRUE, facet by compound name
compounds	numeric vector of compound numbers to include. If NULL, include all compounds

Value

ggplot2 object

plot_var_pattern	<i>Plot Relationship Between Concentration and CV/SD</i>
------------------	--

Description

Plot Relationship Between Concentration and CV/SD

Usage

```
plot_var_pattern(df, title = "")
```

Arguments

df	Data frame with columns: stdconc (standardized concentration), cv (coefficient of variation), sdev (standard deviation), Type (e.g., "Estimated", "Observed")
title	Plot title

Author(s)

Omar I. Elashkar

precision_per_vial	<i>Precision per vial</i>
--------------------	---------------------------

Description

Precision per vial

Usage

```
precision_per_vial(peaks_res, suitability = FALSE)
```

Arguments

peaks_res	PeakRes object
suitability	logical. If TRUE, suitability samples are ignored

Value

ggplot2 object

prefilter_precision_data
Filter data

Description

Filter data

Usage

```
prefilter_precision_data(  
  x,  
  type,  
  acc_cutoff = 0.2,  
  dev_cutoff = 0.2,  
  compound_id = NULL  
)  
  
## S4 method for signature 'QuantRes'  
prefilter_precision_data(  
  x,  
  type,  
  acc_cutoff = 0.2,  
  dev_cutoff = 0.2,  
  compound_id = NULL  
)  
  
## S4 method for signature 'data.frame'  
prefilter_precision_data(x, type, acc_cutoff = 0.2, dev_cutoff = 0.2)
```

Arguments

x	Dataframe or QuantRes Object
type	QC, DQC, or Standard
acc_cutoff	Accuracy cutoff. 20% by default
dev_cutoff	Deviation cutoff. 20% by default
compound_id	Compound ID to filter. If NULL, all compounds are considered

Value

Filtered data

Author(s)

Omar I. Elashkar

`quant_app`*Quantification App*

Description

This function creates a shiny app for quantification after peak integration

Usage

```
quant_app()
```

`read_chrom`*Read Chromatogram Files*

Description

This function reads chromatogram files from a directory and returns a data frame with the chromatogram data.

Usage

```
read_chrom(dir, format = "waters_raw", method)
```

Arguments

<code>dir</code>	directory for chromatograms
<code>format</code>	format of the chromatogram files. Options are "waters_raw" and "mzML".
<code>method</code>	LC-MS/MS method ID saved available in the database.

Examples

```
## Not run:  
path <- system.file("extdata", "waters_raw_ex", package="PKbioanalysis")  
main <- read_chrom(path, method = 1)  
  
## End(Not run)
```

read_experiment_results	<i>Read experiment results</i>
-------------------------	--------------------------------

Description

Read experiment results

Usage

```
read_experiment_results(x, drop_prefix = FALSE, vendor = "targetlynx_xml")
```

Arguments

x	path to experiment results. See details
drop_prefix	logical. If TRUE, drop the prefix from the sample name
vendor	vendor name. Currently only "targetlynx_xml" or "targetlynx_csv" are supported.

Details

Currently only targetlynx XML or CSV exported files are supported.

Value

QuantRes object with the results of the experiment.

register_plate	<i>This will save the plate to the database</i>
----------------	---

Description

This will save the plate to the database

Usage

```
register_plate(plate)
```

Arguments

plate	PlateObj object or MultiPlate object
-------	--------------------------------------

Value

PlateObj object or list of PlateObj objects

response_to_conc	<i>Convert response to concentration</i>
------------------	--

Description

Convert response to concentration

Usage

```
response_to_conc(quantres, compound_id, response)
```

Arguments

quantres	QuantRes object
compound_id	character
response	numeric. Must match the response type used in linearity. Either abs_response or rel_response

Value

numeric

reverse_predict	<i>Reverse predict concentration from response</i>
-----------------	--

Description

Reverse predict concentration from response

Usage

```
reverse_predict(fit, newdata, intercept)
```

Arguments

fit	lm object
newdata	vector or data frame with response values
intercept	logical. Whether the model has intercept or not

Value

numeric. Estimated concentration

Author(s)

Omar I. Elashkar

run_summary	<i>Get Summary of an object</i>
-------------	---------------------------------

Description

Get Summary of an object

Usage

```
run_summary(object)

## S3 method for class 'PeakRes'
run_summary(object)
```

Arguments

object	A PeakRes object
--------	------------------

smooth_chrom	<i>Smooth Chromatogram Peaks</i>
--------------	----------------------------------

Description

This function smooths chromatogram peaks using different algorithms.

Usage

```
smooth_chrom(chrom_res, filter = "mean", window = 2, iter = 2)
```

Arguments

chrom_res	ChromRes object
filter	Filter to use. Options are "mean", "median", "savgol", "gaussian"
window	Window size for the filter
iter	Number of iterations. If 0, no smoothing is applied.

study_app

bioanalytic_app

Description

This function creates a shiny app for plate management

Usage

```
study_app()
```

Value

A shiny app. No default return value. Can return a PlateObj if reuse_plate_button is clicked

update_RT

Manually Update Observed RT for either all compounds, all next samples, or single compound and sample

Description

Update RT for either all compounds, all next samples, or single compound and sample

Usage

```
update_RT(
  chrom_res,
  compound_id,
  sample_id = NULL,
  peak_start,
  peak_end,
  manual = FALSE,
  target = "single",
  force = FALSE
)
```

Arguments

chrom_res	ChromRes object
compound_id	Compound ID
sample_id	Sample ID
peak_start	Minimum RT value
peak_end	Maximum RT value
manual	Manual update. Default is FALSE
target	Target of update. Options are "single", "all", "all_next"
force	Force update if previous peak exists. Default is FALSE

Details

Only target = "all" will update the expected RT for all compounds.

Examples

```
## Not run:
update_RT(chrom_res, compound_id = 1, sample_id = 1,
          peak_start = 1, peak_end = 2, target = "single")

## End(Not run)
```

write_injec_seq	<i>Write injection sequence to database</i>
-----------------	---

Description

Write injection sequence to database

Usage

```
write_injec_seq(injec_seq)
```

Arguments

injec_seq InjecListObj object

Value

dataframe

[[,MultiPlate-method	<i>Subsetting method for MultiPlate</i>
----------------------	---

Description

Subsetting method for MultiPlate

Usage

```
## S4 method for signature 'MultiPlate'
x[[i, j, ...]]
```

Arguments

x	MultiPlate object
i	index
j	index
...	additional arguments

Value

PlateObj object

Index

`[[,MultiPlate-method, 39`

`add_blank, 3`
`add_cs_curve, 4`
`add_DB, 5`
`add_DQC, 5`
`add_QC, 6`
`add_samples, 7`
`add_samples_db, 8`
`add_samples_db2, 8`
`add_suitability, 9`
`area_report.PeakRes, 10`

`calc_var_summary, 10`
`check_chrom_cmpds, 11`
`chrom_app, 12`
`combine_injec_lists, 12`
`combine_plates, 13`
`config_suitability, 13`
`create_new_study, 14`
`cv, 14`

`download_sample_list, 15`

`estim_dil_limit, 15`
`estim_lloq, 16`
`export_integration, 16`
`export_run, 17`
`extract_peak_bounds, 17`

`fill_scheme, 18`
`filter_chrom, 18`
`fit_var, 19`
`formatted_print, 20`

`generate_96, 20`
`get_compound_ID, 21`
`get_sample_ID, 21`
`get_sample_names, 22`

`has_default_RT, 22`

`install_py_dep, 23`
`integrate, 23`
`is_integrated, 24`
`is_integrated,ChromRes-method
(is_integrated), 24`
`is_integrated,ChromResBase-method
(is_integrated), 24`
`is_smoothed, 24`

`length,MultiPlate-method, 25`

`make_calibration_study, 25`
`make_metabolic_study, 26`

`plate_metadata, 27`
`plate_tree, 27`
`plot.PlateObj, 28`
`plot_chrom, 29`
`plot_peak_areas.PeakRes, 30`
`plot_RT.ChromRes, 31`
`plot_RT.PeakRes, 31`
`plot_var_pattern, 32`
`precision_per_vial, 32`
`prefilter_precision_data, 33`
`prefilter_precision_data,data.frame-method
(prefilter_precision_data), 33`
`prefilter_precision_data,QuantRes-method
(prefilter_precision_data), 33`

`quant_app, 34`

`read_chrom, 34`
`read_experiment_results, 35`
`register_plate, 35`
`response_to_conc, 36`
`reverse_predict, 36`
`run_summary, 37`

`smooth_chrom, 37`
`study_app, 38`

update_RT, [38](#)

write_injec_seq, [39](#)